

# PAULDRACH medical GmbH

# SESA



## Absaugpumpe - Suction pump

## 185-300-10

### Applikation:

- Absaugen von Sekreten aus dem oberen und unteren Verdauungstrakt

### Eigenschaften:

- Leistungsstarke Absaugpumpe, 30 l/min
- Äußerst laufruhig
- Glas-Zylinder Kolbensystem
- Pumpengehäuse mit Tragegriff
- Vakuumregler, Manometer
- Hydrophober Bakterienfilter
- Verbindungsschlauch
- 2 Liter Sekretflasche und
- Sekretflaschendeckel mit Griff
- Überlaufsicherung
- Maße: 345 x 245 x 282 mm, 7,4 kg
- Netzspannung 230 V, Schutzklasse 1



### Application:

- Secretion aspirator for the upper and lower intestinal tract.

### Features:

- High suction capacity pump, 30l/min
- Especially low noise level
- Glass-cylinder piston pump system
- Casing with carrying handle
- Vacuum control valve and gauge
- Hydrophobic bacterial filter
- 2 litres suction bottle
- Suction bottle top lid
- Overflow valve
- Dimensions: 345 x 245 x 282 mm,
- Weight: 7.4 kg, 230 Volts

**Absaugpumpe „Senator“**  
komplett mit 2 Liter Behälter,  
Vakuumregler, Manometer,  
Überlaufsicherung,  
Bakterienfilter

**185-300-10**

**Suction pump „Senator“**  
*Comes with 2 liters suction  
bottle, control valve, gauge,  
overflow valve, bacterial filter*

**Hydrophober  
Bakterienfilter**

**185-301-01**

**Hydrophobic bacterial  
filter**

**2 Liter Sekretflasche**

**185-301-02**

*2 liter suction bottle*

**Flaschendeckel mit Griff**

**185-301-03**

*Top cover with handle*

28.08.2004

Distributors stamp

**PAULDRACH medical GmbH**

Porschestr. 22 D- 30827 Garbsen

☎ Inland: 05131-496515 Fax 92014

☎ Intern'l.: +49-5131-496516 Fax 496519

Internet: [www.pauldrach-medical.de](http://www.pauldrach-medical.de)

e-mail: [info@pauldrach.de](mailto:info@pauldrach.de)