

PAULDRACH

medical GMBH

STATE OF THE ART

Reinigungsbürste

13



Produktinformationen

Reinigungsbürste
Bürstenkopf \varnothing 3 mm
mit abgerundeter Spitze
flexibler distaler Spirale
Handgriff + Daumenring
Spirale \varnothing 1,85 x 2300mm

Product information

Cleaning brush, reusable
bristles \varnothing 3 mm
rounded metal tip
flexible distal metal sheath
plastic handle with thumb ring
sheath \varnothing 1.85 x 2300 mm

170-105-30

Reinigungsbürste 3 mm
Metallspirale 180 cm
cleaning brush 3 mm
metal sheath 180 cm

170-108-30

Reinigungsbürste 3 mm
Metallspirale 230 cm
cleaning brush 3 mm
metal sheath 230 cm