

PAULDRACH

medical GMBH

STATE OF THE ART

1 x Biopsie-Zangen

21

1X



Produktinformationen

1 x Biopsie Zange
mit / ohne Dorn, \varnothing 2,3 mm
Metallspirale 2,3 mm
Längen: 180 cm + 230 cm
steril, einzeln verpackt
Mindestbestellmenge: 10 Stück

Product information

Disposable biopsy forceps
with / without needle \varnothing 2.3 mm
metal sheath 2.3 mm
lengths: 180 cm + 230 cm
sterile, single wrapped
minimum order: box of 10

145-915-01

Zange mit Dorn, 180 cm
forceps with needle

145-915-02

Zange ohne Dorn, 180 cm
forceps without needle

155-919-01

Zange mit Dorn, 230 cm
forceps with needle

155-919-02

Zange ohne Dorn, 230 cm
forceps without needle