

# PAULDRACH


## medical GMBH

STATE OF THE ART

### 1 x Polypektomieschlinge

22

1X



#### Produktinformationen

1 x Polypektomieschlinge  
oval, drehbar,  $\varnothing$  25 mm  
Teflonschlauch 2,3 mm  
Längen: 180 cm, 230 cm  
Standard HF-Anschluss  
Längenstopper  
Mindestbestellmenge: 10 Stück

#### Product information

Disposable polypectomy snare  
oval loop, rotatable,  $\varnothing$  25 mm  
teflon sheath 2.3 mm  
lengths: 180 cm, 230 cm  
standard hf-connection  
with stopper  
minimum order: box of 10

147-945-15 180 cm  
ovale Schlinge  $\varnothing$  15 mm  
oval snare  $\varnothing$  15 mm

147-945-25 180 cm  
ovale Schlinge  $\varnothing$  25 mm  
oval snare  $\varnothing$  25 mm

157-949-15 230 cm  
ovale Schlinge  $\varnothing$  15 mm  
oval snare  $\varnothing$  15 mm

157-949-25 230 cm  
ovale Schlinge  $\varnothing$  25 mm  
oval snare  $\varnothing$  25 mm